



PRIONICS AG

WAGISTRASSE 27A
8952 SCHLIEREN-ZURICH
SWITZERLAND

PHONE +41 44 200 20 00
FAX +41 44 200 20 10

INFO@PRIONICS.COM
WWW.PRIONICS.COM

Established Salmonellosis detection. PrioCHECK® SALMONELLA Ab porcine



RELIABLE, SIMPLE, CONVENIENT AND ECONOMICAL TEST FOR SMALL
AND LARGE SCALE SCREENING OF SALMONELLA INFECTIONS IN PIGS.


WWW.PRIONICS.COM

PrioCHECK® Salmonella Ab porcine

Rely on PrioCHECK®



PrioCHECK® Salmonella Ab porcine



PrioCHECK® Salmonella Ab

Prod. No. 7610480 50 plate kit (4500 tests)
Prod. No. 7610460 5 plate kit (450 tests)

SALMONELLOSIS

Salmonellosis is one of the most important zoonotic diseases caused by the organisms of the two species of *Salmonella* (*Salmonella enterica* and *Salmonella bongori*). The presence of Salmonellosis has been identified in countries throughout the world, but appears to be most prevalent in areas of intensive animal husbandry, especially occurring in pigs and calves and some types of poultry. The disease can affect all domestic animals; however young animals, pregnant and lactating animals are the most susceptible to the disease. Many animals, especially pigs and poultry, may also be infected but show no clinical symptoms. Such animals may be important in relation to the spread of infection between flocks and herds and also as a source of food contamination and human infections. In order to control and minimize the risk of contamination, control of the (pork) supply chain is necessary.

Prionics offers an effective tool for the diagnosis of Salmonella in pigs. The PrioCHECK® Salmonella Ab porcine is a reliable, convenient and economical test for small and large scale screening of Salmonella infections in pigs.

WHAT ARE THE BENEFITS OF THE PrioCHECK® SALMONELLA Ab porcine?

Reliable and proven test

The PrioCHECK® Salmonella Ab porcine is the commercial version of the Danish Salmonella mix-ELISA as developed by the Danish Veterinary Institute (DVI). The Danish test has been successfully and reliably applied in the Danish Salmonella control program and is considered to be the gold standard for Salmonella testing. The PrioCHECK® Salmonella Ab porcine is the only test that uses the same critical components (antigen and reference serum).

The PrioCHECK® Salmonella Ab porcine has proven its reliability in Danish control programs and is therefore the best choice in Salmonella control programs..

Convenient

The PrioCHECK® Salmonella Ab porcine follows a simple, 4-step protocol that can easily be automated and connected to a laboratory management system, allowing automated tracking from sample to result. In order to follow long-term variations in handling of the test, a validation control is included on each plate, as was recommended by the OIE.

Its simple and short protocol make the PrioCHECK® Salmonella Ab porcine fully automatable.

Versatile

This versatile test is validated for serum, plasma and meat juice samples and detects antibodies against all serotypes with public health significance (the serogroups B, C1 and D). The PrioCHECK® Salmonella Ab porcine is an economical test for small and large scale screening of Salmonella infections in pigs and is available in small (5 plate) or bulk (50 plate) package sizes.

The PrioCHECK® Salmonella Ab porcine contributes to food safety by detecting all significant Salmonella types significant for human health.

CONTACT INFORMATION

PRIONICS LELYSTAD B.V.
PLATINASTRAAT 33
8211 AR LELYSTAD
THE NETHERLANDS

Tel. +31 320 71 40 00
Fax +31 320 71 40 29
INFO@PRIONICS.NL
WWW.PRIONICS.NL