

One product for every case forensiX



EASY TO USE

- **Sample and swab sealed and dried in the tube** – no folding of a cardboard container required
- **No re-training of personnel required** – similar tube design across the expansive forensiX product line
- **Ensured pairing of swab and tube** – identical ID-Number on both the swab cap and the tube

COST SAVING

- **Avoid costs associated with contaminated swabs** – all parts are DNA free *
- **Reduced invalid/negative PCR results** – fast drying leads to improved sample quality
- **Lower transportation costs** – refrigeration or rapid transport not required

SOPHISTICATED

- **No degradation of DNA after long term storage**
- **Can be used in automated and non-automated labs** – 2D-barcode and ID-number provided on tamper proof label
- **Independent of environmental condition** – drying not influenced by humidity



WWW.FORENSIX.CH

forensiX ONE FOR EVERY CASE



EVIDENCE COLLECTION


When collecting evidence at a crime scene, it is important to secure any usable evidence and observing the chain of custody until completion of analysis.

Superior fast drying of samples in forensiX evidence collection tubes prevents and inhibits growth of bacteria and fungi ensuring protection of valuable DNA samples. Very rapid drying also enables preservation of as much DNA as possible.



The evidence collection kits contain tubes and sterile water to increase the usability even further.

Tubes for evidence collection are available in two lengths: 100mm and 150mm.



BUCCAL SWAB

Build on the same platform as the evidence collection tube, forensiX has created a tube for reference samples. The tube has its own unique color coded label.

Prevention of accidental mixing-up of samples is very important. The tamper proof seal is clearly pairing the swab with the tube, as both, swab and tube have the same unique ID-Number. This number can also be read electronically as a 2D-barcode.



SEXUAL ASSAULT CARE KITS

Collecting and preserving for extended periods all possible evidence in a sexual assault case is critical. The «Hessen Kit», which has been developed in cooperation with the LKA Hessen ¹, is easy to use, reliably secures the evidence and preserves it for years. The kit employs the long self-drying tubes with its own color coded labels.

¹ Hessen State Criminal Office

* DNA-free means: Any DNA, that might contaminate the swab or tube during the manufacturing process, is degraded to such small fragments that no human profile can be generated with PCR amplification methods.

Evidence Collection Tube

Art. Nr. 9022030 150 mm

Art. Nr. 9022035 100 mm

Evidence Collection Kit 150 mm

Art. Nr. 9023130 (5x) 150 mm

Art. Nr. 9023030 (1x) 150 mm

Evidence Collection Kit 100 mm

Art. Nr. 9023135 (5x) 100 mm

Art. Nr. 9023035 (1x) 100 mm

Buccal Swab Tube

Art. Nr. 9012030 150 mm

Hessen Kit

Art. Nr. 9032010

CONTACT INFORMATION

PRIONICS AG
WAGISTRASSE 27A
CH-8952 SCHLIEREN-ZURICH
SWITZERLAND

PHONE +41 44 200 20 00
FAX +41 44 200 20 10
INFO@PRIONICS.COM
WWW.PRIONICS.COM