


DNA collection products for crime scene investigations forensiX




The forensiX product line provides simple to use and highly reliable solutions for collecting DNA samples in criminal investigations. forensiX DNA sampling products are one of the best solutions on the market today because collection, transport, processing, storage and archiving of sample swabs can be carried out at room temperature and without pre-drying.



WWW.SWISSFORENSIX.COM

forensiX

SIMPLE TO USE AND RELIABLE



DNA SAMPLING IN CRIMINAL INVESTIGATIONS

One of the most important tools used for criminal investigative procedures today is DNA analysis and identification. DNA evidence is traceable in extremely small quantities and professional preservation of biological clues obtained from DNA analysis can be critical to resolving criminal cases. Prionics forensiX sampling products conform to the highest quality standards required for the preservation of DNA evidence and are used by many police departments and forensic institutes worldwide.

Prionics acquired the forensiX product lines in July 2008 and now offers a comprehensive range of products for obtaining DNA samples in criminal investigations.




WHAT ARE THE BENEFITS OF THE forensiX PRODUCTS?

forensiX DNA sampling products have been successfully tested in practice and are simple to use. Distinctive to the forensiX product is the drying of cotton swabs in a dedicated storage carton. DNA degradation, one of the biggest risks in forensic sampling, is avoided through the drying process. DNA samples placed in forensiX containers do not require cooling. Transport, storage and archiving are all possible at room temperature. Each storage box is secured with three tamper-proof seals, which are labeled distinctively, allowing an uninterrupted chain of evidence.



Prionics forensiX kits include all sampling materials and detailed instructions for taking DNA swabs. Instructions for use are available in several languages. Please contact us at www.prionics.com or call +41 44 2002000 to find out if we offer the kits in your preferred language.



The **forensiX evidence collection** kit for collecting, packaging and preserving biological evidence for DNA analysis contains 5 sample collection boxes carrying 10 sterile cotton swabs. Five sterile water dispensers are also included in the kit for collection of dry DNA samples.

The **forensiX buccal swab** kit contains 1 buccal swab collection box with 2 sterile cotton swabs to collect samples for DNA analyses. With one kit, sampling of an individual for both DNA analysis and storage in a DNA database can be performed.

The new **forensiX DNA sample storage tubes** can be used either for the collection of criminal evidence samples or for buccal samples. The tubes contain a desiccant that ensures complete and rapid drying of the cotton swab after securing of evidence. The tubes, which were developed together with Dr. Harald Schneider of the state Office of Criminal Investigation of Hessen in Wiesbaden (Germany), fully comply with all known requirements for DNA sampling devices.

forensiX evidence collection kit

Prod. No. 90210..

forensiX buccal swab kit

Prod. No. 90110..

forensiX DNA sample storage tubes for evidence collection

Prod. No. 90220..

forensiX DNA sample storage tubes for buccal swabs

Prod. No. 90120..

CONTACT INFORMATION

PRIONICS AG
WAGISTRASSE 27A
CH-8952 SCHLIEREN-ZURICH
SWITZERLAND

PHONE +41 44 200 20 00
FAX +41 44 200 20 10
INFO@PRIONICS.COM
WWW.PRIONICS.COM